

NEU Proteomics Capstone Module: Case Studies in Quantitative Proteomics

Jacob D. Jaffe, Lev Litichevski, Jarrett Egerston,
Mike MacCoss, Meena Choi, Olga Vitek
w/ special guest: Ruedi Aebersold

Goals for this course

- Understand the concepts and value of a proteomic sentinel assay
- Understand the considerations in moving from discovery proteomics data to targeted proteomics
 - Practical: Selection of probes for a targeted proteomics sentinel assay
- Understand considerations for treatment and normalization of research grade targeted proteomics assays
 - Practical: Construct your own data processing pipeline
- Introduction to next-gen MS methods to expand targeted panels
- Understand the conceptual differences between using markers and profiles of activity
- Introduction to the concept of connectivity and exploration of the dataset in this space

Key Differences from Previous Modules

- This is the first time we have attempted this course
 - Expect some bumps in the road!
- Hackathons, not exercises
 - Put into practice all you have learned so far
 - Team-based work
 - There are no right or wrong answers
 - Key learning comes from sharing of experiences

Course Structure

	Wednesday (May 10)	Thursday (May 11)	Friday (May 12)
8:00-9:00		Finish Hackathon 1	
9:00-9:30		Compare, Contrast, and Discuss Results from Hackthon 1, Compare to P100 Set	Expansion of P100 via Next-Gen MS (DIA)
9:30-10:00			Mini Hackathon: Hack 1 probe observability
10:00-10:30			
10:30-11:00		refreshments	
11:00-11:30		The Targeted P100 Assay in Practice and Considerations for Reproducible Research	The Statistician's View: What's right and what's wrong with the P100 Methodologies
11:30-12:00			
12:00-12:30		Prep for Hackathon 2 - Data intro, code review	
12:30-13:00	Registration	Lunch	
13:00-13:30			
13:30-14:00	Keynote Address: Ruedi Aebersold	Hackathon 2 - Build your own data processing pipeline	Using Sentinel Profiles for Drug Characterization: Clustering and connectivity
14:00-14:30			
14:30-15:00			
15:00-15:30	refreshments		
15:30-16:00	Intro to Sentinel Assay Concept, Development of P100 Discovery Data Set	Keynote: Mike MacCoss - Proteomics Data Ecosystems	Closing Address: Ruedi Aebersold
16:00-16:30			
16:30-17:00	Prep for Hackathon 1 - Data intro, code review		
17:00-17:30	Begin Hackathon 1 - Probe Selection	Compare, Contrast, and Discuss Results from Hackthon 2	Wrap-up
17:30-18:00			

Everything you need:

bit.ly/PCCSEData

The screenshot shows the PanoramaWeb interface for the LINCS PCCSE Overview. The top navigation bar includes the LINCS logo and the title "LINCS PCCSE Overview". A sidebar on the left contains a "LINCS" menu with a tree view of the data structure. The main content area displays a list of items under the "LINCS" heading, including "P100", "GCP", "LINCS_Phase_I", "P100 Template", "GCP Template", "Overview Information", "DIAdev", "AutoQCdev", "Abelin et al - Supplemental Skyline Documents", "Similarity and Connectivity Matrices", and "NEU Capstone Course 2017". A "SITORY" section on the right is partially visible, showing text about "naling and I" and "tions (drug modificatio". At the bottom, a timeline labeled "Time Resolved" shows three stages: "Acute stimulation (small molecule)", "P100", and "RNA". A yellow lightning bolt icon is positioned below the "Acute stimulation" label.

Panorama PanoramaWeb

LINCS LINCS PCCSE Overview

LINCS

- LINCS**
 - P100
 - GCP
 - LINCS_Phase_I
 - P100 Template
 - GCP Template
 - Overview Information
 - DIAdev
 - AutoQCdev
 - Abelin et al - Supplemental Skyline Documents
 - Similarity and Connectivity Matrices
 - NEU Capstone Course 2017

SITORY

naling and I
tions (drug modificatio

Acute stimulation (small molecule) **P100** **RNA**

Time Resolved

Lecture 1: Targeted Proteomics, the Sentinel Assay Concept, and Initial Development of the P100 Assay

Jacob D. Jaffe

Targeted vs. Discovery Proteomics




Shotgun Proteomics

- Complex sample is interrogated stochastically
- Get protein/peptide IDs
- Irreproducible repeats

Targeted Proteomics

- Focus efforts on relatively fewer analytes
- Verify ID, get quantification
- Consistency


Targeted MS: Experimental Implementation



- 100s of analytes in a single experiment
- **Synthetic peptide internal standards** – quantification and proof of ID
- ***When you want to guarantee you measure it each and every time!***

Why do targeted studies?

- Instrumentation inherently undersamples population
 - Speed
 - Sensitivity
 - Dynamic Range



Data Shapes in Proteomics

1. Long and skinny

Peptides/proteins		1	2
AK1	AK1	AK1	AK1
AK2	AK2	AK2	AK2
AK3	AK3	AK3	AK3
AK4	AK4	AK4	AK4
AK5	AK5	AK5	AK5
AK6	AK6	AK6	AK6
AK7	AK7	AK7	AK7
AK8	AK8	AK8	AK8
AK9	AK9	AK9	AK9
AK10	AK10	AK10	AK10
AK11	AK11	AK11	AK11
AK12	AK12	AK12	AK12
AK13	AK13	AK13	AK13
AK14	AK14	AK14	AK14
AK15	AK15	AK15	AK15
AK16	AK16	AK16	AK16
AK17	AK17	AK17	AK17
AK18	AK18	AK18	AK18
AK19	AK19	AK19	AK19
AK20	AK20	AK20	AK20
AK21	AK21	AK21	AK21
AK22	AK22	AK22	AK22
AK23	AK23	AK23	AK23
AK24	AK24	AK24	AK24
AK25	AK25	AK25	AK25
AK26	AK26	AK26	AK26
AK27	AK27	AK27	AK27
AK28	AK28	AK28	AK28
AK29	AK29	AK29	AK29
AK30	AK30	AK30	AK30
AK31	AK31	AK31	AK31
AK32	AK32	AK32	AK32
AK33	AK33	AK33	AK33
AK34	AK34	AK34	AK34
AK35	AK35	AK35	AK35
AK36	AK36	AK36	AK36
AK37	AK37	AK37	AK37
AK38	AK38	AK38	AK38
AK39	AK39	AK39	AK39
AK40	AK40	AK40	AK40
AK41	AK41	AK41	AK41
AK42	AK42	AK42	AK42
AK43	AK43	AK43	AK43
AK44	AK44	AK44	AK44
AK45	AK45	AK45	AK45
AK46	AK46	AK46	AK46
AK47	AK47	AK47	AK47
AK48	AK48	AK48	AK48
AK49	AK49	AK49	AK49
AK50	AK50	AK50	AK50
AK51	AK51	AK51	AK51
AK52	AK52	AK52	AK52
AK53	AK53	AK53	AK53
AK54	AK54	AK54	AK54
AK55	AK55	AK55	AK55
AK56	AK56	AK56	AK56
AK57	AK57	AK57	AK57
AK58	AK58	AK58	AK58
AK59	AK59	AK59	AK59
AK60	AK60	AK60	AK60
AK61	AK61	AK61	AK61
AK62	AK62	AK62	AK62
AK63	AK63	AK63	AK63
AK64	AK64	AK64	AK64
AK65	AK65	AK65	AK65
AK66	AK66	AK66	AK66
AK67	AK67	AK67	AK67
AK68	AK68	AK68	AK68
AK69	AK69	AK69	AK69
AK70	AK70	AK70	AK70
AK71	AK71	AK71	AK71
AK72	AK72	AK72	AK72
AK73	AK73	AK73	AK73
AK74	AK74	AK74	AK74
AK75	AK75	AK75	AK75
AK76	AK76	AK76	AK76
AK77	AK77	AK77	AK77
AK78	AK78	AK78	AK78
AK79	AK79	AK79	AK79
AK80	AK80	AK80	AK80
AK81	AK81	AK81	AK81
AK82	AK82	AK82	AK82
AK83	AK83	AK83	AK83
AK84	AK84	AK84	AK84
AK85	AK85	AK85	AK85
AK86	AK86	AK86	AK86
AK87	AK87	AK87	AK87
AK88	AK88	AK88	AK88
AK89	AK89	AK89	AK89
AK90	AK90	AK90	AK90
AK91	AK91	AK91	AK91
AK92	AK92	AK92	AK92
AK93	AK93	AK93	AK93
AK94	AK94	AK94	AK94
AK95	AK95	AK95	AK95
AK96	AK96	AK96	AK96
AK97	AK97	AK97	AK97
AK98	AK98	AK98	AK98
AK99	AK99	AK99	AK99
AK100	AK100	AK100	AK100


2. Short and wide

[illegible]


3. Hole-y

[illegible]



It's hard to be hole-y



Instrument Jockeys



Data Geniuses



The Sentinel Assay Concept

A Sentinel Assay for Yeast Biology

A sentinel protein assay for simultaneously quantifying cellular processes

Martin Soste¹, Rita Hrabakova^{1,2}, Stefanie Wanka³,
Andre Melnik¹, Paul Boersema¹, Alessio Maiolica⁴,
Timon Wernas¹, Marco Tognetti¹, Christian von Mering³
& Paola Picotti¹

Nature Methods 11:1045-8 (2014)

ously, in high throughput. The approach we describe here is based on the concept of 'sentinel' proteins that report on the activation state of a specific biological process and are selected on the basis of literature evidence or computational prediction (Fig. 1). Our

- Two targeted proteomics assays:
 - Protein
 - Phosphopeptide
- Characterizes a range of biochemical responses to stimulation
- Targets nominated by literature, prior knowledge

A Sentinel Assay for Yeast Biology






Fig. 1 from Soste et al., *Nature Methods* 11:1045-8 (2014)


Examples of Specific Sentinel Behaviors




Selected from Fig. 3 from Soste et al., *Nature Methods* 11:1045-8 (2014)

Examples of Sentinel Reporting Profiles

Protein Level Reporters



Phospho Level Reporters



Selected from Fig. 2 from Soste et al., *Nature Methods* 11:1045-8 (2014)

Considerations for Translating to a Human Sentinel Assay

- Human proteome/phosphoproteome is bigger
- Pathway architecture in human cells is more complex
- Probably less direct knowledge as to what the sentinels would be
- Proteomic accessibility of desired sentinel peptides
- Diversity of biology and cell types coupled with proteome dynamic range might make life difficult

“P100” – A different kind of sentinel assay


Reduced Representation Phosphosignaling

Why phosphoproteomics?

- Classic “pathways” are commonly represented as signaling cascades of phosphoproteomic events
- Phosphoproteomic events are the “first responders” of cells and occur on a range of time scales from minutes to days
- Phosphoproteins are not easily quantified using nucleic acid based techniques
- Kinase inhibitors are a major class of emerging therapeutics

Reduced Representation Signaling Set: P100

- Phosphosignaling is likely coordinated
 - Kinase-substrate relationships tend to be 1-to-many
 - *A priori* expectation of coordinate regulation of sites



- Don't need to monitor every phosphosite all the time?
 - Amenable to reductionist approaches?

Final goals of P100 assay development


- Drug responsiveness
- Universal (-ish)
- Ser/Thr focused
- Compact
- Operates at scale
- Longitudinally comparable

Major Problem #1: No easily usable prior data


- Many piecemeal, *ad hoc* drug studies
 - 1 or 2 drugs
 - Many focus on pTyr signaling
 - Wide range of biological models
 - Inconsistent data reporting formats
 - General lack of systematic studies
- We decided to make our own
 - Discovery proteomics
 - Quantitation by SILAC

Selection of Compounds for Systematic Data Generation


- Start with gene expression profiles induced by drug treatment
 - Multiple cell types
 - Extract kinase/phosphatase genes



Generation of foundational data




Select representative phosphopeptide probes from coherent clusters



- Chosen based on:
- Modulates kinase/phosphatase expression
 - Concentration from CMAP
 - Have matching gene expression data
 - Several same class examples
 - PosCon = staurosporine


Jinal Patel, Xiaodong Lu

SILAC and Stable Isotopes: Let cells do the work



- Feed cells growth medium with “heavy” amino acids
- Proteins incorporate these amino acids
- All proteins are slightly “heavier” but otherwise biochemically equivalent
- Choice of Lys and Arg typical because of trypsin

SILAC-based discovery profiling





- Cells pre-labeled by SILAC
- Treatments: 6 hours, concentrations from gene expression study
- IMAC-based phosphopeptide enrichment
 - Primary yield: Ser/Thr
- Mass spec provides both ID (MS2) and quantification (MS1)

Considerations for successful assay configuration

- Are there a set of phosphosites that we can (almost) always observe?
- Do these phosphosites behave in interesting ways in response to perturbations?
- Do we gain anything over gene expression measurements?
- Can we pick a (small) set of highly informative phosphopeptides?
- Could profiling these peptides reproduce the observations of a deeper experiment?
- Is the assay generally useful?


Experimental Implementation



More than 7.5×10^8 cells grown for experiments!
True bio-replicates

Computational Fractionation Acquisition Strategy


- All experiments performed on a Q-Exactive mass spectrometer
 - 2 hour, “Top 12” shotgun proteomics experiments
 - Each sample analyzed 3 times:



- Primary interpretation and quantification in MaxQuant

P100 Source Dataset

Jinal Patel, Xiaodong Lu




- Large systematic set of perturbations with phosphoproteomics readout
 - Over 10,000 phosphosites observed
- Over 1,200 sites present in >75% of all experiments

Will it work?



- ✓ Are there a set of phosphosites that we can (almost) always observe?
- Do these phosphosites behave in interesting ways in response to perturbations?
- Do we gain anything over gene expression measurements?
- Can we pick a (small) set of highly informative phosphopeptides?
- Could profiling these peptides reproduce the observations of a deeper experiment?
- Is the assay generally useful?

Signals span lineages




- All 6 Taxol treatments cluster – regardless of lineage


Signals group related compounds and replicates




Digoxigenin




Digitoxigenin



Digoxin



Lanatoside C




Anisomycin

Will it work?

- ✓ Are there a set of phosphosites that we can (almost) always observe?
- ✓ Do these phosphosites behave in interesting ways in response to perturbations?
- Do we gain anything over gene expression measurements?
- Can we pick a (small) set of highly informative phosphopeptides?
- Could profiling these peptides reproduce the observations of a deeper experiment?
- Is the assay generally useful?

Phospho-state adds value over gene expression

A. Information angle between phosphosignaling vs. gene expression





Will it work?

- ✓ Are there a set of phosphosites that we can (almost) always observe?
- ✓ Do these phosphosites behave in interesting ways in response to perturbations?
- ✓ Do we gain anything over gene expression measurements?
- Can we pick a (small) set of highly informative phosphopeptides?
- Could profiling these peptides reproduce the observations of a deeper experiment?
- Is the assay generally useful?

Determining the Right Number of Groups to Monitor


- Each group should contain phosphopeptides with similar patterns
- PCA and HC methods suggest ~50 clusters are optimal
- Further refinement by hierarchical clustering leads to 55 clusters
- We can choose 2 probes per cluster to build redundancy into assay



Principles for picking group representatives


1. # sites (1 is best)
2. Missed cuts (0 is best)
3. C and M count (0 is best)
4. Observability (156 is best)
5. Localization (1 is best)
6. Length (shorter is better)

Dimensionality reduction: >1000 to 55 with 2x probe redundancy




Probeset recalls replicates significantly better than random

Random Probes



P100 probes



- Looks like it is going to work.

Will it work?

- ✓ Are there a set of phosphosites that we can (almost) always observe?
- ✓ Do these phosphosites behave in interesting ways in response to perturbations?
- ✓ Do we gain anything over gene expression measurements?
- ✓ Can we pick a (small) set of highly informative phosphopeptides?
- Could profiling these peptides reproduce the observations of a deeper experiment?
- Is the assay generally useful?

To be continued...

- But first, a hackathon